

ILLUMINATING INNOVATION: VERSATILE FLUORESCENT MICROPARTICLES FOR DYNAMIC DETECTION

Mason's innovative fluorescent microparticles offer a groundbreaking approach to sensor technology and molecular detection. Utilizing a one-pot Reversible Addition-Fragmentation Chain Transfer (RAFT) polymerization method, researchers have created versatile, reversibly quenchable, fluorescent hydrogel microparticles. These particles are temperature-resistant and offer robust, long-lasting fluorescence, making them suitable for a range of applications, from environmental sensing to membrane coverage validation.

Key Features

- **One-Pot Synthesis:** Simple and efficient RAFT polymerization process to create molecularly imprinted, fluorescent microparticles
- **Reversible Quenching:** Fluorescent microparticles can be quenched and restored, allowing for dynamic detection and sensor applications
- **Temperature Stability:** Particle size remains stable across a temperature range of 25°C to 50°C, with minimal changes in fluorescence even under long-term storage conditions
- **Versatility:** These microparticles can be tailored for various applications, including analyte capture, environmental monitoring, and bioimaging
- **Multiple Fluorophores:** Easily adaptable to different fluorophores without the need for reengineering, offering flexibility in detection strategies
- **Cost-Effective:** Offers a scalable and efficient solution for producing fluorescent particles without significant modifications to the synthesis process

These unique fluorescent microparticles provide a customizable and adaptable platform for numerous scientific and industrial applications, with reliable fluorescence for precise detection and analysis.



For More Information contact:

George Mason University, Office of Technology Transfer
703-993-8933 ott@gmu.edu <https://ott.gmu.edu/>